



IASTM Tools

IASTM-Tools
Outils IASTM
Herramientas IASTM
Strumenti IASTM



Mark Approved



20+

Years of
TIG
Performance



Welcome to T.I.G International.

T.I.G International. is a knowledge and research based organization, which was established in 2002 in Sialkot (Pakistan), a city which is globally famous for high quality Surgical Instruments and IASTM Tools manufacturing.

Our main focus is on the design and development of instruments for Minimally Invasive Surgery and our current product range involves Reusable and Single Use Surgical Instruments , IASTM Tools and training and education solutions for surgical Skill Labs. New products are regularly added to our ever-increasing products portfolio, a proof our firm commitment on research and development.

Our dedicated technical staff works closely with doctors to understand their needs and then convert their ideas into practical products with high index of safety and performance. Our company puts alot of emphasis on customer satisfaction. We achieve this through an excellent products quality, quick delivery and most importantly, an efficient after-sales service.

Our Vision:

To be recognized as Pakistan's No.1 Surgical Solution provider.

Our Mission:

To improve patient safety be promoting surgical education through providing cheaper and beeta surgical solutions.

We look forward to being of service!

To promote surgical education.

T.I.G International Teams.

CONTENTS

| Description | Page No |
|-----------------------------|---------|
| IASTM TOOL MUSTACHE..... | 001 |
| IASTM TOOL KISS..... | 002 |
| IASTM TOOL WAVE..... | 003 |
| IASTM TOOL ARTIST..... | 004 |
| IASTM TOOL SET LARGE..... | 005 |
| IASTM TOOL RAZOR..... | 006 |
| IASTM TOOL SET SMALL..... | 007 |
| IASTM HANDLE BAR LARGE..... | 008 |
| IASTM HANDLE BAR SMALL..... | 009 |
| IASTM TOOL DROP..... | 010 |
| IASTM TOOL GRIP..... | 011 |
| IASTM TOOL WHEELER..... | 012 |

Note: We have mentioned the standard sizes and weight in our catalogue. If you have required the your own size and weight, so we will make the products as per your sizes and weight.

Reusable Tools



TIT-001

The T.I.G International® IASTM Tool – Mustache is a wide applicable, light to hold and our best selling IASTM Tool. The “Mustache” has a flat finish on the bottom. The top has both round and upwardly slanted edges. The instrument also has a pointy side to reach deeper fascial layers. A popular choice for general use. Travel pouch included.

Material: AISI 304 (1.4301) high-quality surgical stainless steel

Weight: 0.33 lbs

Dimensions: 11 × 3 × 1.2 in

Reusable Tools



TIT-002


The T.I.G International® IASTM Tool – Kiss has rounded edges on both sides and a firm hold. Extremely suitable for stroking in both directions. Also well-suited for other techniques. Typical for this tool is its long, round, convex curve. Travel pouch included.

Material: AISI 304 (1.4301) high-quality surgical stainless steel

Weight: 300 grams

Dimensions: 5mm x 9cm x 17cm

Reusable Tools



TIT-003

The T.I.G International® IASTM TOOL IASTM Tool – Wave was designed to act as an all-in-one IASTM tool. One tool allowing for many different application techniques. Travel pouch included.

Reusable Tools



TIT-004


The T.I.G International® IASTM Tool – Artist is a robust IASTM Tool with a straight finish at the bottom and an inward-faced flat edge at the top. Ideal for scanning tissue. This Tool has different shapes, such as a concave side, a convex side and jagged corners. This makes it widely applicable. Travel pouch included.

Material: AISI 304 (1.4301) high-quality surgical stainless steel

Weight: 240 grams

Dimensions: 4mm x 11cm x 14cm

Reusable Tools



The T.I.G International® IASTM Tool Set size Large consists of five large tools for treating all parts of the body and larger muscle groups. Travel pouch included.

The set size L consists of five large tools for treating all parts of the body and larger muscle groups.

Material: AISI 304 (1.4301) high-quality surgical stainless steel, Total weight of 1310 grams.

Type 1: (9mm x 1.9cm x 44.6cm)
Weight: 562 grams

Type 2: (6mm x 6cm x 13.6cm)
Weight: 186 grams

Type 3: (6mm x 4.2cm x 20.7cm)
Weight: 228 grams

Type 4: (6mm x 3.2cm x 21.3cm)
Weight: 193 grams

Type 5: (5mm x 2.3cm x 15.2cm)
Weight: 78 grams

Reusable Tools



TIT-006


The T.I.G International® IASTM Tool – Razor is a comfortable to hold IASTM Tool. The tool has different round and flat finishes across the entire tool on one side. It is widely applicable for treatment, scanning and stroking. Travel pouch included.

Material: AISI 304 (1.4301) high-quality surgical stainless steel

Weight: 325 grams

Dimensions: 5mm x 3.8cm x 19.2cm

Reusable Tools



TIT-007

The IASTM Tool Set size Small consists of five small, narrow, fine tools, particularly for treating smaller parts of the body. Travel pouch included. Total weight of 315 grams.

Small Set of 5 Instruments

Material: AISI 304 (1.4301) high-quality surgical stainless steel

Type 1: (5mm x 1cm x 31cm)
Weight: 144 grams

Type 2: (3mm x 1.8cm x 15.5cm)
Weight: 59 grams

Type 3: (3mm x 1.8cm x 12.6cm)
Weight: 46 grams

Type 4: (3mm x 3cm x 11cm)
Weight: 35 grams

Type 5: (3mm x 1.4cm x 14.5cm)
Weight: 41 grams

Reusable Tools



TIT-008

High quality T.I.G International® IASTM tool for treatment of bigger surfaces. Two handles and curved in the middle for deep tissue scraping. Travel pouch included.

Material: AISI 304 (1.4301) high-quality surgical stainless steel

Weight: 0.539 kg

Dimensions: 44.6 × 1.9 × 0.9 cm

Weight: 562 grams

Reusable Tools



TIT-009

T.I.G International® IASTM tool with comfortable handles and rounded edges. Therefore, this tool is particularly useful for longer scrapes. Travel pouch included.

Material: AISI 304 (1.4301) high-quality surgical stainless steel

Weight: 0.153 kg

Dimensions: 31 × 1 × 0.5 cm

Weight: 144 grams

Reusable Tools



TIT-010

The T.I.G International® IASTM Tool – Drop is a light IASTM Tool with a flat round finishing. Because the “Drop” is only 3 mm thick, it is ideal for stroking in both directions. Also suitable for other fascial techniques. Travel pouch included.

Material: AISI 304 (1.4301) high-quality surgical stainless steel

Weight: 138 grams

Dimensions: 3mm x 6.5cm x 16cm

Reusable Tools



TIT-011

Due to its various convex and concave finished 'double beveled' treatment surfaces, the T.I.G International® Grip IASTM tool is suitable for treating a wide range of areas. This includes fingers as well as larger surfaces like the legs. Travel pouch included.

Material: AISI 304 (1.4301) high-quality surgical stainless steel

Weight: 0.9 lbs

Dimensions: 9 × 6.2 × 1.2 inch

Reusable Tools



TIT-012

The T.I.G International® IASTM Tool – Wheeler is a generally applicable flat IASTM Tool (4 mm) with various corners and sides (identical on both sides). The “Wheeler” also has a button to reach deeper fascial layers. Suitable for various fascial techniques (scanning, stroking and treatment). Travel pouch included.

Material: AISI 304 (1.4301) high-quality surgical stainless steel

Weight: 0.51 lbs

Dimensions: 9 × 6.1 × 1 inch

CONTACT US:



T.I.G International
3/199 Street No.1,
Puran Nagar Paris Road,
Sialkot-51310, Pakistan.

Tel: +92 52 3559212
Fax: +92 52 3559212
Hotline: +92 300 6166818
WhatsApp: +92 333 8785863
Email: info@tig-int.net / tigint786@hotmail.com
URL: www.tig-int.net

